


SELECTION OF COUPLING FOCAL LENGTH FOR AN ELLIPTICAL BEAM

DETERMINATION OF THE EFFECTIVE COUPLING DIAMETER

In order to find the best coupling focal length in case of an elliptical beam use the effective beam diameter \varnothing_{eff} which is calculated from the small and the large diameters $\varnothing_{||}$ and \varnothing_{\perp} of the collimated elliptical laser beam:

$$\varnothing_{eff} = \sqrt{\varnothing_{||} \cdot \varnothing_{\perp}}$$

Definition of the effective beam diameter \varnothing_{eff} .



SELECTION OF FOCAL LENGTH

After determining the effective beam diameter \varnothing_{eff} . The best focal length can be determined. This is done by replacing the diameter \varnothing_{beam} with \varnothing_{eff} in the corresponding [formula](#).



5AN-2.5

MEASURES TO INCREASE THE COUPLING EFFICIENCY FOR ELLIPTICAL BEAMS

Coupling efficiency can be increased by converting the elliptical laser beam to a circular beam using [anamorphic beam-shaping optics](#).

This is a printout of the page

<https://sukhamburg.com/support/technotes/fiberoptics/coupling/couplingsm/focallengthelliptical.html> from 5/22/2026

CONTACT

For more information please contact:

Schäfter + Kirchhoff GmbH

Kieler Str. 212

22525 Hamburg

Germany

Tel: +49 40 85 39 97-0

Fax: +49 40 85 39 97-79

info@sukhamburg.com

www.sukhamburg.com

LEGAL NOTICE

Copyright 2020 Schäfter+Kirchhoff GmbH. All rights reserved.

Text, image, graphic, sound, video and animation files and their arrangement on Schäfter+Kirchhoff GmbH webpages are protected by copyright and other protective laws. The content may not be copied for commercial use or reproduced, modified or used on other websites. [\[more\]](#)